

Resource Summary Report

Generated by [NIF](#) on Jul 31, 2026

Peroxidase-AffiniPure Rabbit Anti-Mouse IgG + IgM (H+L) (min X Hu Sr Prot)

RRID:AB_2340069

Type: Antibody

Proper Citation

Jackson ImmunoResearch Labs Cat# 315-035-048, RRID:AB_2340069

Antibody Information

URL: http://antibodyregistry.org/AB_2340069

Proper Citation: Jackson ImmunoResearch Labs Cat# 315-035-048, RRID:AB_2340069

Target Antigen: Mouse IgG + IgM (H+L)

Clonality: unknown

Comments: Originating manufacturer of this product

Antibody Name: Peroxidase-AffiniPure Rabbit Anti-Mouse IgG + IgM (H+L) (min X Hu Sr Prot)

Description: This unknown targets Mouse IgG + IgM (H+L)

Antibody ID: AB_2340069

Vendor: Jackson ImmunoResearch Labs

Catalog Number: 315-035-048

Record Creation Time: 20250819T235114+0000

Record Last Update: 20250820T063718+0000

Ratings and Alerts

No rating or validation information has been found for Peroxidase-AffiniPure Rabbit Anti-Mouse IgG + IgM (H+L) (min X Hu Sr Prot).

No alerts have been found for Peroxidase-AffiniPure Rabbit Anti-Mouse IgG + IgM (H+L) (min X Hu Sr Prot).

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 10 mentions in open access literature.

Listed below are recent publications. The full list is available at [NIF](#) .

Gupta D, et al. (2025) COP9 signalosome and PRMT5 methylosome complexes are essential regulators of Lis1-dynein-based transport. *Cell reports*, 45(1), 116736.

Gupta D, et al. (2025) A LisH-domain protein interaction map reveals a Lis1-ARIH2-dynein regulatory axis. *iScience*, 28(11), 113912.

Reshi HA, et al. (2024) EYA protein complex is required for Wntless retrograde trafficking from endosomes to Golgi. *Developmental cell*, 59(18), 2443.

Bülow S, et al. (2024) Bactericidal/permeability-increasing protein instructs dendritic cells to elicit Th22 cell response. *Cell reports*, 43(3), 113929.

Xiao M, et al. (2024) Smad4 sequestered in SFPQ condensates prevents TGF- β tumor-suppressive signaling. *Developmental cell*, 59(1), 48.

Iriki T, et al. (2023) Senescent cells form nuclear foci that contain the 26S proteasome. *Cell reports*, 42(8), 112880.

Holzinger JM, et al. (2023) Scorpionfish BPI is highly active against multiple drug-resistant *Pseudomonas aeruginosa* isolates from people with cystic fibrosis. *eLife*, 12.

Dai L, et al. (2021) Structural insight into BRCA1-BARD1 complex recruitment to damaged chromatin. *Molecular cell*, 81(13), 2765.

Koyano F, et al. (2019) Parkin-mediated ubiquitylation redistributes MITOL/March5 from mitochondria to peroxisomes. *EMBO reports*, 20(12), e47728.

Uesato T, et al. (2018) Enhanced Expression of the Key Mitosis Regulator Cyclin B1 Is Mediated by PDZ-Binding Kinase in Islets of Pregnant Mice. *Journal of the Endocrine Society*, 2(3), 207.