

Resource Summary Report

Generated by [NIF](#) on Sep 11, 2026

Thermo Fisher: Q Exactive Hybrid Quadrupole-Orbitrap Mass Spectrometer

RRID:SCR_020565

Type: Tool

Proper Citation

Thermo Fisher: Q Exactive Hybrid Quadrupole-Orbitrap Mass Spectrometer
(RRID:SCR_020565)

Resource Information

URL:

<https://www.thermofisher.com/order/catalog/product/IQLAAEGAAPFALGMAZR#/IQLAAEGAAPFALGMAZR>

Proper Citation: Thermo Fisher: Q Exactive Hybrid Quadrupole-Orbitrap Mass Spectrometer (RRID:SCR_020565)

Description: Identify, quantify and confirm more compounds rapidly and with confidence using the Thermo Scientific Q Exactive Hybrid Quadrupole-Orbitrap Mass Spectrometer. This benchtop LC-MS/MS system combines quadrupole precursor ion selection with high-resolution, accurate-mass (HRAM) Orbitrap detection to deliver performance and versatility. The Q Exactive Mass Spectrometer is equally useful for untargeted or targeted screening and a broad range of qualitative and quantitative applications in drug discovery, proteomics, environmental and food safety, clinical research and forensic toxicology.

Resource Type: instrument resource

Keywords: Thermo, Thermo LTQ, Thermo Scientific, Thermo Exactive, Thermo Fisher, Thermo Trace, Thermo Q, Thermo Electrom, Mass Spectrometer, Instrument Equipment, USEdit

Availability: Commercially available

Resource Name: Thermo Fisher: Q Exactive Hybrid Quadrupole-Orbitrap Mass Spectrometer

Resource ID: SCR_020565

Alternate IDs: Model_Number_Q Exactive Hybrid

Record Creation Time: 20220129T080351+0000

Record Last Update: 20260905T132920+0000

Ratings and Alerts

No rating or validation information has been found for Thermo Fisher: Q Exactive Hybrid Quadrupole-Orbitrap Mass Spectrometer.

No alerts have been found for Thermo Fisher: Q Exactive Hybrid Quadrupole-Orbitrap Mass Spectrometer.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We found 3 mentions in open access literature.

Listed below are recent publications. The full list is available at [NIF](#) .

K??cher S, et al. (2026) Overexpression of the ERG oncogene in prostate cancer identifies candidates for PARP inhibitor-based radiosensitization. The Journal of clinical investigation, 136(6).

Ravazza D, et al. (2026) Quantitative Mass Spectrometry-Based Biodistribution of Monoclonal Antibodies: An Alternative to Radio-Biodistribution. Advanced science (Weinheim, Baden-Wurttemberg, Germany), 13(38), e20265.

Preza M, et al. (2022) Global analysis of neuropeptides in cestodes identifies Attachin, a SIFamide homolog, as a stimulant of parasite motility and attachment. Journal of neurochemistry, 162(6), 467.